
F&P Optiflow™ Junior 2

Next generation care

Nasal Cannula INSTRUCTIONS FOR USE

OJR410

OJR412

OJR414

OJR416


OJR418

Manufacturer Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348, Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.co.nz Web: www.fphcare.com
Australia (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 Fax: +61 3 9871 4998 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554
China 代理人/售后服务机构: 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号
 电话: +86 20 32053486 传真: +86 20 32052132 **Finland** Tel: +358 (0)405 406618 Fax: +46 (0)8 36 6310 **France** Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 **Germany** Tel: +49 7181 98599 0 Fax: +49 7181 98599 66 **India** Tel: +91 80 4284 4000 Fax: +91 80 4123 6044 **Irish Republic** Tel: 1800 409 011 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea** Tel: +82 2 6205 6900 Fax: +82 2 6309 6901 **Northern Ireland** Tel: 0800 132 189 **Russia** Tel and Fax: +7 495 782 21 50
Spain Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **Taiwan** Tel: +886 2 8751 1739 Fax: +886 2 8751 5625 **Turkey** Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, Alinteri Bulvarı 1161/1 Sokak No. 12-14, P.O. Box 06371 Ostim, Ankara, Turkey Tel: +90 312 354 34 12 Fax: +90 312 354 31 01
UK Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: +44 1628 626 136 Fax: +44 1628 626 146 **USA/Canada** Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

SIZING INFORMATION

1. Recommended nare occlusion of 50% should be used to size cannula.
2. Patient weight should only be used as a guide.

		WEIGHT GUIDE (KG)											
		500g	1	1.5	2	2.5	3	3.5	4	SPARES			
	OJR410 XS	█		█		█						F&P Wigglepads™ 2 WJR110	
	OJR412 S	█		█		█		█		█			

		WEIGHT GUIDE (KG)															
		1	2	3	4	5	6	7	8	9	10	15	20	25	30	SPARES	
	OJR414 M	█		█		█		█		█						F&P Wigglepads™ 2 WJR112	
	OJR416 L	█		█		█		█		█		█		█			
	OJR418 XL	█		█		█		█		█		█		█			

3. Expected to fit patient May fit patient










TECHNICAL SPECIFICATIONS

FLOW RATES (L/MIN)

Platform	MR850	AIRVO™ 2
Circuit	RT330 RT331	900PT531 900PT561*
OJR410 XS	0.5-8	N/A
OJR412 S	0.5-9	
OJR414 M	0.5-10	2-20
OJR416 L	0.5-23	
OJR418 XL	0.5-25	

4. Flow rates above describe technical capability of the product when used at sea level. Ensure clinical judgement is used when prescribing flow rates.
5. Maximum flow rates achieved with the RT331 are ventilator dependant. Refer to the manufacturer's instructions for correct usage.

*Not available for sale in the United States.

I Apply Cannula**II Remove Cannula****III Replace F&P Wigglepads™ 2****F&P Optiflow™ Junior 2**

Next generation care

F&P Optiflow™ Junior 2 nasal cannula**Intended Use**

The Fisher & Paykel Healthcare Optiflow™ Junior 2 nasal cannula is a single use nasal cannula intended for use with a nasal high flow therapy system to deliver heated and humidified nasal high flow therapy to spontaneously breathing patients who require breathing support.

This product is designed for use in hospital environments and must be prescribed by a physician.

The intended pediatric subpopulations targeted for use of the F&P Optiflow™ Junior 2 nasal cannula range includes:

- Neonates, birth up to 1 month of age
- Infants, 1 month up to 2 years of age
- Children, 2 years up to 12 years of age

General Warnings

- This product is only designed and verified for use with equipment, accessories and spare parts approved by F&P. Unauthorized equipment, accessories or spare parts which are used with this product may impair performance of this product or compromise safety (including potentially causing serious patient harm).
- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in the event of an interruption to gas flow) may result in serious harm or death.
- If using supplemental oxygen, keep ignition sources away from the patient.

GENERAL CAUTIONS

- Regularly monitor the patient to ensure skin integrity and that the skin underneath the cannula remains dry. A barrier film may be used between the cannula and the patient's upper lip to prevent irritation.
- DO NOT soak, sterilize, or reuse this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers. Secretions on the cannula and the prongs can be removed by gently wiping with a damp cloth.
- Reuse may result in the transmission of infectious substances, interruption to treatment, serious harm or death.
- DO NOT stretch the cannula on application; this may cause increased pressure to the patient's skin. If necessary, the cannula may be repositioned.
- DO NOT use in an MRI or similar magnetic-based scanner as the tubing contains stainless steel.
- Tubing may pose a risk of strangulation or airway restriction.
- DO NOT use if the product or its packaging has been tampered with.
- Ensure the patient does not lie on the tubing as this may apply pressure to the patient's ears or face.
- Product to only be used with medical grade gas supplies. The gas supply used with this device may unexpectedly fail to deliver oxygen or flow.
- Failure to apply and use this product within the directions, transport, storage and operating conditions specified in the labeling and user instructions may impair performance of this product or compromise safety (including potentially causing serious patient harm).

Optiflow, Wigglepads, and the Seahorse, Crab, Starfish, Octopus and Turtle images shown in this user instruction are trade marks of Fisher & Paykel Healthcare Limited.

© 2017 Fisher & Paykel Healthcare Limited.

For patent information, see www.fphcare.com/ip

California residents please be advised of the following, pursuant to Proposition 65:

WARNING

This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm.

For more information, please visit: www.fphcare.com/prop65

I Apply Cannula

- 1 Select appropriate cannula size; recommended nare occlusion is approximately 50%.
 - Prepare the patient's skin according to hospital protocol.
 - Connect the system to the gas source and ensure there is gas flow through the prongs.
- 2 Remove the first backing tabs from the F&P Wigglepads™ 2 and avoid touching the adhesive.
- 3 Insert the cannula into the nares. Ensure the cannula bridge rests close to the nose without touching the septum. Do not stretch the cannula during application. Stick the Wigglepads 2 to the patient's cheeks.
- 4 Remove the second backing tabs and stick the Wigglepads 2 onto the cheeks.

WARNINGS

- DO NOT allow the prongs to seal in the nares. Occlusion may result in septal damage or barotrauma.
- When using the RT331 breathing circuit, a ventilator with Nasal High Flow mode and appropriate pressure limits must be used. Failure to do so may compromise therapy, leading to serious patient harm (e.g. hypoxia or barotrauma).

CAUTIONS

- DO NOT place the Wigglepads 2 on the patient's eyes, ears or injured skin.
- Ensure the cannula is placed directly on the Wigglepads 2. Direct skin contact caused by cannula misalignment may result in skin breakdown.

II Remove Cannula

Place fingertip on the outside edge of the Wigglepads 2 and gently peel the cannula away from the Wigglepads 2. Starting from the outside, peel towards the nose.

III Replace F&P Wigglepads™ 2

- Lift the edge of the Wigglepads 2. Use a damp cloth to wipe the patient's skin and the underside of the Wigglepads 2 while gently peeling away from the patient's face.
 - 1 Adhere the replacement Wigglepads 2 to the underside, remove the first backing tabs and stick the Wigglepads 2 onto the patient's cheeks.
 - 2 Remove the second backing tabs and stick onto the patient's cheeks.

Checks During Operation

- Regular monitoring of the patient is necessary to ensure a slight gap between the cannula and septum is maintained, as well as the correct placement of the prongs in the nares. Reposition the cannula on the Wigglepads 2 if required.
- Regularly check skin integrity to prevent excessive pressure to the upper lip.
- To prevent nare occlusion, clean secretions from the cannula and the patient's face as required.
- Check the cannula remains secure. Replace the Wigglepads 2 if required.
- Ensure all connections are secure during use. Check the cannula is undamaged and that the flow path is maintained. Under excessive load, the cannula may disconnect to prevent forces being transferred to the patient.

CAUTIONS

- DO NOT wrap, insulate, stretch or crush the tubing as this may impair the performance of this product or compromise safety (including potentially causing patient harm).
- Monitor condensate to prevent occlusion or build-up of fluid. Drain away from the patient as required.

Contraindications

This therapy must not be used where continuous positive airway pressure (CPAP) is contraindicated. This includes:

- Non-spontaneous breathing.
- Injury, congenital abnormalities and anatomical malformations where bi-nasal prongs or physical conditions where positive airway pressure is contraindicated, including but not limited to; pneumothorax, pneumocephalus, cerebrospinal fluid leak and hypotension.
- Injury/trauma/severe deformity that might be exacerbated by use of nasal prongs or mask.

Technical Specifications

This product is single use.

This product is not made with natural rubber latex.

This product does not contain PVC or Phthalates (DEHP, DBP, BBP).

This product is intended to be used for a maximum of 7 days.

CAUTION

Using this product beyond 7 days may impair performance of this product or compromise safety (including potentially causing serious patient harm).

OPERATING CONDITIONS

- Ambient temperature range: 18 to 26 °C.

NOTE

This product may be used in an incubator.

**APPROVED COMPATIBLE EQUIPMENT/
ACCESSORIES**

- MR850 Humidifier in invasive mode with approved breathing circuits, accessories and chamber kits.
- AIRVO™ 2 Humidifier with approved breathing circuits, accessories and chamber kits.
- OPT014 Oxygen Tubing.

Refer to the table for a complete list of approved breathing circuits and accessories.

OPERATING FLOW RATES

The operating flow rates for each F&P Optiflow™ Junior 2 cannula size are dependent on the breathing circuit and humidifier in use.

Refer to page 2 for the operating flow rates for each cannula.

- OPT014 Oxygen Tubing flow rates 0.1-2 L/min for all cannula sizes.

SIZING NOTES

1. Recommended nare occlusion of 50% should be used to size cannula.
2. Patient weight should only be used as a guide.
3.  Expected to fit patient.  May fit patient.

TECHNICAL SPECIFICATION NOTES

4. Flow rates above describe technical capability of the product when used at sea level. Ensure clinical judgement is used when prescribing flow rates.
5. Maximum flow rates achieved with the RT331 are ventilator dependant. Refer to the manufacturer's instructions for correct usage.

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