Aerogen[®]

Pioneering Aerosol Drug Delivery

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11

Références

- 1 Ari A, Atalay OT, Harwood R, Sheard MM, Aljamhan EA and Fink JB. Influence of nebulizer type, position, and bias flow on aerosol drug delivery in simulated pediatric and adult lung models during mechanical ventilation. *Respiratory care*. 2010;55:845-51.
- 2 Berlinski A and Willis JR. Albuterol delivery by 4 different nebulizers placed in 4 different positions in a pediatric ventilator in vitro model. *Respiratory care*. 2013;58:1124-33.
- 3 Fang TP, Lin HL, Chiu SH, Wang SH, DiBlasi RM, Tsai YH and Fink JB. Aerosol Delivery Using Jet Nebulizer and Vibrating Mesh Nebulizer During High Frequency Oscillatory Ventilation: An In Vitro Comparison. *Journal of aerosol medicine and pulmonary drug delivery*. 2016.
- 4 Abdelrahim ME, Plant P and Chrystyn H. In-vitro characterisation of the nebulised dose during non-invasive ventilation. *The Journal of pharmacy and pharmacology*. 2010;62:966-72.
- 5 Galindo-Filho VC, Ramos ME, Rattes CS, Barbosa AK, Brandao DC, Brandao SC, Fink JB and de Andrade AD. Radioaerosol Pulmonary Deposition Using Mesh and Jet Nebulizers During Noninvasive Ventilation in Healthy Subjects. *Respiratory care*. 2015;60:1238-46.
- 6 McPeck M. Improved Aerosol Drug Delivery with an Electronic Mesh Nebulizer during Non-invasive Ventilation *poster de l'AARC*. 2012
- 7 Ari A, Harwood R, Sheard M, Dailey P and Fink JB. In vitro comparison of heliox and oxygen in aerosol delivery using pediatric high flow nasal cannula. *Pediatric pulmonology*. 2011;46:795-801.
- 8 Reminiac F, Vecellio L, Heuze-Vourc'h N, Petitcollin A, Respaud R, Cabrera M, Le Pennec D, Diot P and Ehrmann S. Aerosol Therapy in Adults Receiving High Flow Nasal Cannula Oxygen Therapy. *Journal of aerosol medicine and pulmonary drug delivery*. 2015.
- 9 Alcoforado L, de Melo Barcelar J, Castor Galindo V, Cristina S, Brandão S, Fink J, B. and Dornelas de Andrade A. Analysis of Deposition Radioaerosol Nebulizers Membrane in Healthy Subjects. *Poster de présentation de l'ISAM* 2015.
- 10 Ari A, de Andrade A, Sheard M, AlHamad B and Fink JB. Performance Comparisons of Jet and Mesh Nebulizers Using Different Interfaces in Simulated Spontaneously Breathing Adults and Children. *Journal of aerosol medicine and pulmonary drug delivery*. 2014.
- 11 Dugernier J, Reyçhler G, Depoortere V, Roeseler J, Michotte JB, Laterre PF, Jamar F and Hesse M. Tomoscintigraphic comparison of lung deposition with a vibrating-mesh and a jet nebulizer. *Poster de conférence de l'ERS* 2015.
- 12 Hickin S, Mac Loughlin R, Sweeney L, Tatham A and Gidwani S. Comparison of mesh nebuliser versus jet nebuliser in simulated adults with chronic obstructive pulmonary disease. *Poster à la College of Emergency Medicine Clinical Excellence Conference*. 2014

- 13 Dubus JC, Vecellio L, De Monte M, Fink JB, Grimbert D, Montharu J, Valat C, Behan N and Diot P. Aerosol deposition in neonatal ventilation. *Pediatric research*. 2005;58:10-4.
- 14 Ramsi MA, Henry M, Milla CE and Cornfield DN. Inhaled beta2-Agonist Therapy Increases Functional Residual Capacity in Mechanically Ventilated Children With Respiratory Failure. *Pediatric critical care medicine : a journal of the Society of Critical Care Medicine and the World Federation of Pediatric Intensive and Critical Care Societies*. 2015;16:e189-93.
- 15 Blake G, Yaklic J and Cobb J. Transition from ipratropium/albuterol inhaler to nebulizer on quality and cost savings in ventilated patients. *Poster de présentation à l'American Society of Health-System Pharmacists*. 2013.
- 16 McDaniel C, Glynn G, Gudowski S, Pezzano T and Weibel S. Conversion of pressurized Metered-Dose Inhaler to Vibrating Mesh Nebulizer Administered Medications. *Poster de présentation à la Society of Critical Care Medicine*. 2013.
- 17 Streepy KS, Dawson AM, Grigonis AM, Hammerman SI and Snyder LK. Conversion from metered dose inhalers to a vibrating mesh nebulizer in long term acute care hospitals: cost effectiveness and respiratory staff perception *Poster de présentation à l'American Society of Health-System Pharmacists* 2013.
- 18 Loborec SM, Johnson SE and Keating EA. Financial effect of converting ipratropium-albuterol therapy from inhalers to nebulizer treatments at an academic health system. *American journal of health-system pharmacy : AJHP : journal officiel de l'American Society of Health-System Pharmacists*. 2016;73:121-5.
- 19 Ari A, de Andrade AD, Sheard M, AlHamad B and Fink JB. Performance Comparisons of Jet and Mesh Nebulizers Using Different Interfaces in Simulated Spontaneously Breathing Adults and Children. *Journal of aerosol medicine and pulmonary drug delivery*. 2015;28:281-9.
- 20 Reminiac F, Vecellio L, Heuze-Vourc'h N, Petitcollin A, Respaud R, Cabrera M, Pennec DL, Diot P and Ehrmann S. Aerosol Therapy in Adults Receiving High Flow Nasal Cannula Oxygen Therapy. *Journal of aerosol medicine and pulmonary drug delivery*. 2016;29:134-41.
- 21 Ari A and Fink JB. Guidelines for aerosol devices in infants, children and adults: which to choose, why and how to achieve effective aerosol therapy. *Expert review of respiratory medicine*. 2011;5:561-72.
- 22 Ari A, Areabi H and Fink JB. Evaluation of aerosol generator devices at 3 locations in humidified and non-humidified circuits during adult mechanical ventilation. *Respiratory care*. 2010;55:837-44.
- 23 Sidler-Moix AL, Di Paolo ER, Dolci U, Berger-Gryllaki M, Cotting J and Pannatier A. *Physicochemical aspects and efficiency of albuterol nebulization: comparison of three aerosol types in an in vitro pediatric model*. *Respiratory care*. 2015;60:38-46.
- 24 Ari A, Fink JB and Dhand R. Inhalation Therapy in Patients Receiving Mechanical Ventilation: An Update. *Journal of aerosol medicine and pulmonary drug delivery*. 2012.
- 25 Diot P, Morra L and Smaldone GC. Albuterol delivery in a model of mechanical ventilation. Comparison of metered-dose inhaler and nebulizer efficiency. *American journal of respiratory and critical care medicine*. 1995;152:1391-4.

- 26 Everard ML, Devadason SG, Summers QA and Le Souef PN. Factors affecting total and “respirable” dose delivered by a salbutamol metered dose inhaler. *Thorax*. 1995;50:746-9.
- 27 Emberger J. Cost Reduction using Aeroneb Solo in a Medical ICU Ventilator Population. . *Respiratory Care, Christiana Care Health System, Newark ; abstracts de l'OPEN FORUM 2011*. 2011.
- 28 White CC, Crotwell DN, Shen S, Salyer J, Yung D, Zheng J and DiBlasi RM. Bronchodilator delivery during simulated pediatric noninvasive ventilation. *Respiratory care*. 2013;58:1459-66.
- 29 AlQuaimi M, Fink JB, Harwood R, Sheard M, Bryant L and Ari A. The effect of nebulizer type and mask design on aerosol delivery during noninvasive positive pressure ventilation of an adult lung model. *Poster de l'AARC 2012*.
- 30 Alalwan M, Ari A, Fink J, Harwood R, Bryant L and Sheard M. Delivery of albuterol by pressurized metered-dose inhaler and jet nebulizer via mask with high flow nasal cannula in place reduces aerosol delivery. *Abstract de l'Open Forum de l'AARC*. 2012.
- 31 Bhashyam AR, Wolf MT, Marcinkowski AL, Saville A, Thomas K, Carcillo JA and Corcoran TE. Aerosol delivery through nasal cannulas: an in vitro study. *Journal of aerosol medicine and pulmonary drug delivery*. 2008;21:181-8.
- 32 MacLoughlin R, Power P, Wolny M and Duffy C. Evaluation of vibrating mesh nebulizer performance during nasal high flow therapy. *Poster de présentation à l'International Society of Aerosols in Medicine 2013*.
- 33 Alcoforado L, Ari A, De Melo Barcelar J, Brandao SS, Fink JB and Dornelas De Andrade A. Comparison of Aerosol Deposition with Heated and Unheated High Flow Nasal Cannula (HFNC) in Healthy Adults. *Poster de présentation à l'ATS*. 2016.
- 34 Dailey P, Tina T, Santos J and Gurung P. Performance Improvement Plan for Pediatric Patients in Respiratory Distress: Clinical Experience. *Respiratory Therapy*. 2015;10:27-29.
- 35 Newsome T, Thayer T and Lieberman S. Pediatric Patient with Severe Asthma Treated with a Vibrating Mesh with Valved Adapter. *Respiratory Therapy Case Study*. 2015;10:28-29.
- 36 Thayer T, Hodgens E and Lieberman S. A Case Study: Use of Vibrating Mesh with a Valved Adapter in a Pediatric Patient with a Severe Asthma Exacerbation. *Respiratory Therapy Case Study*. 2016;11:52.
- 37 Dunne R, Shortt SA and Dailey PA. Comparison of bronchodilator administration with vibrating mesh with adapter and jet nebulizer in the emergency department *AARC Open Forum Abstract*. 2016.
- 38 Golden LR, Desimone HA, Yeroshalmi F, Pranevicius M and Saraghi M. Severe intraoperative bronchospasm treated with a vibrating-mesh nebulizer. *Anesthesia progress*. 2012;59:123-6.

Aerogen[®]

Pioneering Aerosol Drug Delivery

Aerogen Ltd.

Galway Business Park,
Dangan, Galway,
Irlande.

INTL. +353 91 540 400

ÉTATS-UNIS (886) 423-7643

www.aerogen.com

PM349