



Home

About Masimo

Signal Extraction  
Pulse Oximetryrainbow® Pulse  
CO-Oximetry

Root

rainbow Acoustic  
Monitoring®Capnography &  
Multigas Monitoring

Brain Monitoring

Patient SafetyNet  
Remote MonitoringOEM Products  
& Solutions

Clinical Evidence

News &amp; Media

Investors


Global Services

Contact

Capnography & Multigas Monitoring Technology | ISA Capnography | ISA OR+ | EMMA Capnometer | EMMA Capnograph | Nomoline

# EMMA™ Mainstream Capnometer

Real-time Capnometry at Your Fingertips



Proud Sponsor of the  
**Patient Safety**  
Science & Technology  
Summit



**\$1 MILLION**  
Masimo SET®  
Performance  
Guarantee

- > EMMA requires minimal warm-up time — just 15 seconds to measure end-tidal carbon dioxide (EtCO<sub>2</sub>) and respiration rate (RR)
- > Small, portable capnometer — lightweight design fits in the palm of your hand for mobility and convenience during short-term EtCO<sub>2</sub> monitoring of adult, pediatric, and infant patients
- > Designed to fit easily onto a breathing circuit
- > Flexible use at multiple environments including pre-hospital, emergency medicine, operating room, and intensive care
- > Rugged design for reliable operation in challenging environments
- > Easy to maintain — no routine calibration required

## Features

- > Simple, easy-to-use interface for quick setup and one-touch programming
- > Audible and visual alarm system for No Breath Detected, No Adapter, Check Adapter, and adjustable high and low EtCO<sub>2</sub> alarm

### One-touch Alarm Limits

One-touch alarm limits and visual alarm status indicator, silence active alarms for two minutes

### Bar Graph

Highly visible bar graph provides continuous feedback on end-tidal carbon dioxide concentration, breathing activity, or alarm status

### Power Button

Warm-up time to full accuracy in 15 seconds

### EtCO<sub>2</sub>

Quantitative EtCO<sub>2</sub> is updated every breath (also available in kPa)

### RR

Respiratory Rate is displayed after two breaths and updated every breath

### Airway Adapter

Comes in adult/pediatric and infant sizes





2014 MASIMO ANNUAL REPORT

thank you



DOWNLOAD NOW

### Kit



> EMMA Kit\*

EMMA (mmHg) PN 9632  
EMMA (kPa) PN 9633

### Accessories



> EMMA Airway Adapter

Adult/Pediatric  
Box of 25  
PN 17448



> EMMA Airway Adapter

Infant  
Box of 10  
PN 17449

\* In order for the EMMA Kit to provide readings, either of the listed airway adapters is required. Kit includes EMMA, pouch, and lanyard.

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

### Resources



**Product Information**  
EMMA Mainstream Capnometer



**Whitepaper**  
Emerging Use of Capnometry in  
Emergency Medicine