


# New Optiflow™+ tracheostomy and mask interfaces

Introducing the OPT970 and OPT980

## OPT970 tracheostomy interface

The new OPT970 tracheostomy interface improves upon the current OPT570 and OPT870 interfaces, featuring a new sputum guard designed with stability and ease of use in mind.



## OPT980 mask interface adapter

The new OPT980 mask interface adapter improves upon the current RT013 interface, for connecting to face and tracheostomy mask interfaces with a standard 22 mm medical taper.



Both new interfaces feature:

### An improved Evaqua™ membrane

Only Fisher & Paykel Healthcare products have Evaqua technology, which reduces the formation of mobile condensate near the patient. Our continuous improvement has led to the newest generation of Evaqua being **stronger** and **more durable**, for clinician confidence when delivering the therapy; and **quieter**, for patient comfort when sleeping. Furthermore, the new Evaqua tubing is **50 mm longer**, providing patients with more freedom of movement.

### A new connector

We have made changes to the connector that fits to the heated breathing tube, making it smaller and less obtrusive, while still ensuring that it provides a secure fit for delivering optimal therapy.

**Please contact your local representative for further information about this or any other Fisher & Paykel Healthcare product or therapy.**