

Adult and Infant Respiratory Care



Our Story

At Fisher & Paykel Healthcare, we fully understand the life-changing impact that can result from delivering inspired and world-leading healthcare solutions. After all, we changed the way the world delivers respiratory humidification to critically ill patients in hospital, and were the first organisation in the world to make a humidifier specifically for use in the treatment of patients with obstructive sleep apnea.

Today, through an unrivalled commitment to innovation; to doing the right thing; and doing what is best for the patient, we continue to inspire world-leading healthcare solutions. Our team of self-motivated, empathetic and committed individuals is our strength. It is this team; our originality; the value we place on relationships; our international perspective; and the utmost care for all patients, customers, suppliers, shareholders, as well as the environment, that will ultimately ensure we will be the healthcare company of choice.

Our Purpose

Improving care and outcomes through inspired and world-leading healthcare solutions.



Welcome



Fisher & Paykel Healthcare is a global leader in design, manufacturing and marketing of medical devices and systems for use in respiratory care and acute care.

Our products and systems are sold in over 120 countries worldwide. We sell our products through direct sales offices operations in most of our major markets, and a network of distributors that sell to hospitals, healthcare providers and other manufacturers of medical devices.

Fisher & Paykel entered the respiratory care market in 1971 with the development of a unique respiratory humidifier system for use in critical care. We now offer a broad range of products and systems for use in respiratory and acute care and in the treatment of OSA.

Contents

AIRVO™ 2 System

Humidifier & accessories	A - 2
--------------------------------	-------

F&P 850™ System B - 1

Humidifier & accessories	B - 2
Humidification chambers	B - 3

F&P 850™ System - Adult

Adult nasal high flow breathing circuits	B - 4
Adult nasal high flow interfaces & accessories	B - 4
Adult noninvasive ventilation breathing circuits	B - 5
Adult face masks & accessories	B - 6
Adult invasive ventilation breathing circuits	B - 7
Adult invasive ventilation interfaces & accessories	B - 8

F&P 850™ System - Infant

Infant resuscitation circuits & kits	B - 9
Infant resuscitator & masks	B - 10
Infant resuscitation accessories	B - 11
Infant invasive ventilation breathing circuits	B - 10
Infant invasive ventilation accessories	B - 11
Infant CPAP therapy breathing circuits	B - 12
Infant CPAP therapy interfaces	B - 13
Infant CPAP therapy accessories	B - 15
Infant nasal high flow breathing circuits	B - 16
Infant nasal high flow interfaces	B - 16
Infant nasal high flow accessories	B - 17

F&P 810™ System

Humidifier & accessories	C - 2
Evatherm™ breathing circuits	C - 2
Humidifier brackets & hardware	D - 1

General humidifier brackets & hardware

Humidifier brackets & hardware	D - 2
Humidifier accessories for reusable systems	E - 1

F&P Reusable Tubing Systems

Reusable humidification chambers	E - 1
Reusable breathing circuit kits	E - 1
Adult reusable breathing circuit (900MR784) spare parts	E - 2
Infant reusable breathing circuit (900MR782) spare parts	E - 3

Disclaimer: Please note that the information in this catalogue (including the product information and images) is summarized and provided for illustrative purposes only. Please refer to the relevant user instructions for more information and confirm details with your local Fisher & Paykel Healthcare representative prior to placing an order. Information is subject to change without notice.

AIRVO™ 2 Humidifier

with Integrated Flow Generator



The AIRVO 2 generates high flows of warmed and humidified respiratory gases, delivered to the patient through a variety of nasal, tracheostomy and mask interfaces.

AIRVO™ 2 System

Product Code	Description	Quantity	
PT101XX	AIRVO™ 2 humidifier with integrated flow generator (XX indicates regional variation)	1 each	
OPT942	Optiflow™+ nasal cannula for AIRVO™/F&P 850 System - Small	20 pack	
OPT944	Optiflow™+ nasal cannula for AIRVO™/F&P 850 System - Medium	20 pack	
OPT946	Optiflow™+ nasal cannula for AIRVO™/F&P 850 System - Large	20 pack	
900PT501	Tube and chamber kit for AIRVO™	10 pack	
OPT842	Optiflow™ nasal cannula for AIRVO™/F&P 850 System - Small	20 pack	
OPT844	Optiflow™ nasal cannula for AIRVO™/F&P 850 System - Medium	20 pack	
OPT846	Optiflow™ nasal cannula for AIRVO™/F&P 850 System - Large	20 pack	
OPT870	Optiflow™ tracheostomy interface for AIRVO™/F&P 850 System	20 pack	
RT013	Mask interface adaptor	20 pack	
900PT531	Tube and chamber kit for AIRVO™ 2 (Junior)	10 pack	
OPT316	Optiflow™ Junior nasal cannula - Infant (Wigglepad™ replacement = OPT012)	20 pack	
OPT318	Optiflow™ Junior nasal cannula - Pediatric (Wigglepad™ replacement = OPT012)	20 pack	

AIRVO™ 2 System

Product Code	Description	Quantity	
900PT422	Oxygen inlet extension kit	1 each	
900PT600	Disinfection kit	1 each	
900PT601	Disinfection filter	2 pack	
900PT602	Cleaning sponge stick	20 pack	
900PT603	Clean storage cover	20 pack	
900PT913	Air filter	2 pack	
900PT405	Pole mounting tray	1 each	
900PT411	UPS mounting kit	1 each	
900PT421	Mobile pole stand	1 each	
900PT426	Accessories basket	1 each	
900PT427	Oxygen cylinder holder	1 each	

F&P 850™ System














The MR850 (F&P 850 System) Heated Humidifier is part of a complete system designed to deliver optimally humidified gas (37°C and 44 mg/L).

Temperature is measured at both the chamber and at the end of the breathing circuit and a dual feedback system incorporates advanced algorithms designed to achieve optimal humidity levels, while minimizing caregiver attention.












The F&P 850 System incorporating the MR850 Heated Humidifier can be used for the full range of Fisher & Paykel Healthcare therapies across the F&P Respiratory Care Continuum™, including invasive and noninvasive ventilation and Optiflow™ therapies.




Humidifier & accessories

Product Code	Description	Quantity	
MR850XXX	MR850 respiratory heated humidifier (XXX indicates regional plug variation)	1 each	
900MR817	MR850 dual heated connection kit	1 each	
900MR805	Heater wire adaptor for dual-heated breathing circuits	1 each	
900MR806	Heater wire adaptor for inspiratory heated breathing circuits	1 each	
900MR863	Temperature/Flow probe adaptor for use with circuits 0.7 m (28") long	1 each	
900MR868	Temperature/Flow probe adaptor for use with circuits 1.1 m (44") long	1 each	
900MR860	Temperature/Flow probe adaptor for use with circuits 1.3 m (52") long	1 each	
900MR869	Temperature/Flow probe adaptor for use with circuits 1.5m (60") long	1 each	
900MR861	Temperature/Flow probe adaptor for use with circuits 1.8m (72") long	1 each	
900MR870	MR850 calibration reference probe	1 each	
900MR208	Reflective shield for temperature probe	20 pack	









Humidification chambers

Product Code	Description	Quantity	
MR290V	Autofeed adult/infant humidification chamber - Vented	40 pack	
MR290VX	Autofeed adult/infant humidification chamber - Vented	10 pack	
MR290HFV	Autofeed adult/infant humidification chamber - High frequency vented	40 pack	
MR210	Manual fill adult humidification chamber - Male/Female ports	40 pack	
MR210X	Manual fill adult humidification chamber - Male/Female ports	10 pack	
MR250	Manual fill adult humidification chamber - Male/Male ports	40 pack	
MR250X	Manual fill adult humidification chamber - Male/Male ports	10 pack	
MR225	Manual fill infant/pediatric humidification chamber	50 pack	
MR225X	Manual fill infant/pediatric humidification chamber	10 pack	
900MR190	Chamber water feed set for MR225/MR250	30 pack	
900MR191X	Chamber water feed set extension for MR290	10 pack	





Adult nasal high flow breathing circuits

Product Code	Description	Quantity	
RT202	Adult inspiratory heated breathing circuit with autofeed chamber	10 pack	




Adult nasal high flow interfaces & accessories

Product Code	Description	Quantity	
OPT942	Optiflow™+ nasal cannula for AIRVO™/F&P 850 System - Small	20 pack	
OPT944	Optiflow™+ nasal cannula for AIRVO™/F&P 850 System - Medium	20 pack	
OPT946	Optiflow™+ nasal cannula for AIRVO™/F&P 850 System - Large	20 pack	
OPT842	Optiflow™ nasal cannula for AIRVO™/F&P 850 System - Small	20 pack	
OPT844	Optiflow™ nasal cannula for AIRVO™/F&P 850 System - Medium	20 pack	
OPT846	Optiflow™ nasal cannula for AIRVO™/F&P 850 System - Large	20 pack	
OPT542	Optiflow™ nasal cannula for F&P 850 System - Small	20 pack	
OPT544	Optiflow™ nasal cannula for F&P 850 System - Medium	20 pack	











Adult nasal high flow interfaces & accessories

Product Code	Description	Quantity	
OPT546	Optiflow™ nasal cannula for F&P 850 System - Large	20 pack	
OPT870	Optiflow™ tracheostomy interface for AIRVO™/F&P 850 System	20 pack	
OPT570	Optiflow™ tracheostomy interface for F&P 850 System	20 pack	
900MR961	Dual flow adaptor - 22 mm Female	20 pack	






Adult noninvasive ventilation breathing circuits

Product Code	Description	Quantity	
RT319	Adult bilevel/CPAP inspiratory heated breathing circuit with autofeed chamber, filter, pressure line & exhalation port	10 pack	
RT301	Adult inspiratory heated breathing circuit with autofeed chamber, Y-piece, CPAP extension, adjustable PEEP valve & pressure relief valve	10 pack	
RT201	Adult inspiratory heated breathing circuit with autofeed chamber, Y-piece & CPAP extension	10 pack	





Adult face masks & accessories

Product Code	Description	Quantity	
RT040L	FreeMotion™ vented full face mask with anti-asphyxiation valve - Large	10 pack	
RT040M	FreeMotion™ vented full face mask with anti-asphyxiation valve - Medium	10 pack	
RT040S	FreeMotion™ vented full face mask with anti-asphyxiation valve - Small	10 pack	
RT040	FreeMotion™ vented full face mask without anti-asphyxiation valve (Australia only)	10 pack	
RT041L	FreeMotion™ non-vented full face mask - Large	10 pack	
RT041M	FreeMotion™ non-vented full face mask - Medium	10 pack	
RT041S	FreeMotion™ non-vented full face mask - Small	10 pack	
RT042	FreeMotion™ vented nasal mask	10 pack	
RT042L	FreeMotion™ vented nasal mask - Large	10 pack	
RT043L	FreeMotion™ non-vented full face mask with anti-asphyxiation valve - Large	10 pack	










Adult face masks & accessories

Product Code	Description	Quantity	
RT043M	FreeMotion™ non-vented full face mask with anti-asphyxiation valve - Medium	10 pack	
RT043S	FreeMotion™ non-vented full face mask with anti-asphyxiation valve - Small	10 pack	
RT016	Inspiratory filter	20 pack	
RT022	Unheated flexible extension	20 pack	
RT023	Patient ventilator weaning kit	20 pack	












Adult invasive ventilation breathing circuits

Product Code	Description	Quantity	
RT385	Adult Evaqua™ 2 dual heated breathing circuit with autofeed chamber, pressure line & filter	10 pack	
RT380	Adult Evaqua™ 2 dual heated breathing circuit with autofeed chamber & filter	10 pack	
RT200	Adult dual heated breathing circuit with autofeed chamber	10 pack	
RT204	Adult dual heated breathing circuit with autofeed chamber & pressure lines	10 pack	








Adult invasive ventilation interfaces & accessories

Product Code	Description	Quantity	
RT019	Inspiratory/Expiratory filter	20 pack	
RT020	Expiratory filter	20 pack	
RT021	DirectConnect™ dual swivel catheter mount	20 pack	
RT023	Patient ventilator weaning kit	20 pack	
RT024	Dry line extension kit	20 pack	
RT028	Circuit hanger	20 pack	
RT032	Water trap	20 pack	
RT913	Inspiratory elbow adaptor - PB7200	10 pack	
RT914	Inspiratory adaptor - PB840	10 pack	










Infant resuscitation circuits & kits

Product Code	Description	Quantity	
900RD110	Infant inspiratory heated resuscitation circuit with T-piece - Humidified	10 pack	
900RD010	Infant T-piece resuscitation circuit with adjustable PEEP valve - Classic	10 pack	
RD1300-10	Infant T-piece resuscitation circuit with adjustable PEEP valve - Ergonomic	10 pack	
900RD014-10	Infant resuscitation kit including T-piece circuit (Classic) & premature mask (S)	10 pack	
900RD015-10	Infant resuscitation kit including T-piece circuit (Classic) & neonatal mask (M)	10 pack	
900RD016-10	Infant resuscitation kit including T-piece circuit (Classic) & infant mask (L)	10 pack	
RD1330-10	Infant resuscitation kit including T-piece circuit (Ergonomic) & micro-premature mask (XS)	10 pack	
RD1340-10	Infant resuscitation kit including T-piece circuit (Ergonomic) & premature mask (S)	10 pack	
RD1350-10	Infant resuscitation kit including T-piece circuit (Ergonomic) & neonatal mask (M)	10 pack	
RD1360-10	Infant resuscitation kit including T-piece circuit (Ergonomic) & infant mask (L)	10 pack	
RD1370-10	Infant resuscitation kit including T-piece circuit (Ergonomic) & pediatric mask (XL)	10 pack	

Infant resuscitator & masks

Product Code	Description	Quantity	
RD900XXX	Neopuff™ infant T-piece resuscitator (XXX indicates fascia language variation)	1 each	
900IW130X	Neopuff™ infant T-piece resuscitator (infant warmer module) (X indicates fascia language variation)	1 each	
RD803-10	Neonatal resuscitation mask - Micro-premature (XS) 35 mm	10 pack	
RD804-10	Neonatal resuscitation mask - Premature (S) 42 mm	10 pack	
RD805-10	Neonatal resuscitation mask - Neonatal (M) 50 mm	10 pack	
RD806-10	Neonatal resuscitation mask - Infant (L) 60 mm	10 pack	
RD807-10	Neonatal resuscitation mask - Pediatric (XL) 72 mm	10 pack	







Infant resuscitation accessories

Product Code	Description	Quantity	
RD020-01	Reusable test lung	1 each	
900RD008	Short gas supply line - 31cm (12.2")	10 pack	
900RD009	Long gas supply line - 210cm (82.7")	10 pack	
900RD101	Gas inlet adaptor	1 each	
900MR148	Reusable T-piece neonatal circuit adaptor - 15 Male/10 Female to 10 Male	1 each	
NC100	NeoWrap™ occlusive wrap	2 pack (NC100-SK) 10 pack (NC100-10) 25 pack (NC100-25)	
900RD301	Side mounting block for Neopuff™	1 each	
900RD302	Wall mounting bracket for Neopuff™	1 each	
RD050-01	Central mounting block kit for Neopuff™	1 each	






Infant invasive ventilation breathing circuits

Product Code	Description	Quantity	
RT265	Infant Evaqua™ 2 dual heated breathing circuit with autofeed chamber & pressure line (flows > 4L/min)	10 pack	
RT266	Infant Evaqua™ 2 dual heated breathing circuit with autofeed chamber & pressure line (flows 0.3 - 4L/min)	10 pack	
RT267	Infant Evaqua™ 2 dual heated breathing circuit with autofeed chamber & pressure line (for SLE2000)	10 pack	
RT268	Infant Evaqua™ 2 dual heated breathing circuit with autofeed chamber & pressure line (for SLE5000)	10 pack	
RT225	Infant inspiratory heated breathing circuit with autofeed chamber & pressure line (flows > 4L/min)	10 pack	
RT226	Infant inspiratory heated breathing circuit with autofeed chamber & pressure line (flows 0.3 - 4L/min)	10 pack	
RT227	Infant inspiratory heated breathing circuit with autofeed chamber & pressure line (for SLE2000/SLE2000 HFO ventilators)	10 pack	
RT228	Infant inspiratory heated breathing circuit with autofeed chamber & pressure line (for SLE4000/SLE5000 ventilators)	10 pack	

Infant invasive ventilation accessories

Product Code	Description	Quantity	
RT010	Inline nebulizer adaptor kit	20 pack	
RT030	Side-stream nebulizer adaptor kit	20 pack	
RT061	Adaptor kit for Dräger (2)	50 pack	
RT062	Luerlock adaptor (2)	50 pack	
RT063	Pressure port adaptor (2)	100 pack	
RT064	Adaptor kit for oxygen/CPAP (2)	50 pack	
RT065	Adaptor kit for Hamilton/EME	10 pack	

Infant CPAP therapy breathing circuits

Product Code	Description	Quantity	
RT324	Infant inspiratory heated breathing circuit with autofeed chamber, pressure line, adaptor kit, expiratory extension & water feed extension (flows > 4L/min)	10 pack	
RT224	Infant inspiratory heated breathing circuit with autofeed chamber, pressure line & adaptor kit (flows > 4L/min)	10 pack	
BC151-10	Infant Bubble CPAP™ system including inspiratory heated breathing circuit with autofeed chamber, pressure manifold & CPAP generator (connects to T-bar interface)	10 pack	
BC161-10	Infant Bubble CPAP™ system including inspiratory heated breathing circuit with autofeed chamber, pressure manifold & CPAP generator (connects to F&P interface)	10 pack	
BC171-10	Infant Bubble CPAP™ system including inspiratory heated breathing circuit with autofeed chamber, pressure manifold with dryline, pressure line & CPAP generator (connects to Flow Driver)	10 pack	






Infant CPAP therapy interfaces

Product Code	Description	Quantity	
BC190-05	FlexiTrunk™ nasal tubing universal fitting (50mm)	5 pack	
BC191-05	FlexiTrunk™ nasal tubing universal fitting (70mm)	5 pack	
BC192-05	FlexiTrunk™ nasal tubing universal fitting (100mm)	5 pack	
BC3020-10	Nasal prong (3.0mm nare diameter / 2.0mm septal width)	10 pack	
BC3520-10	Nasal prong (3.5mm nare diameter/ 2.0mm septal width)	10 pack	
BC4030-10	Nasal prong (4.0mm nare diameter / 3.0mm septal width)	10 pack	
BC4540-10	Nasal prong (4.5mm nare diameter / 4.0mm septal width)	10 pack	
BC5040-10	Nasal prong (5.0mm nare diameter / 4.0mm septal width)	10 pack	
BC5050-10	Nasal prong (5.0mm nare diameter / 5.0mm septal width)	10 pack	
BC5550-10	Nasal prong (5.5mm nare diameter / 5.0mm septal width)	10 pack	
BC5560-10	Nasal prong (5.5mm nare diameter / 6.0mm septal width)	10 pack	
BC6060-10	Nasal prong (6.0mm nare diameter / 6.0mm septal width)	10 pack	
BC6070-10	Nasal prong (6.0mm nare diameter / 7.0mm septal width)	10 pack	
BC6570-10	Nasal prong (6.5mm nare diameter / 7.0mm septal width)	10 pack	



Infant CPAP therapy interfaces

Product Code	Description	Quantity	
BC800-10	Infant nasal mask - Small	10 pack	
BC801-10	Infant nasal mask - Medium	10 pack	
BC802-10	Infant nasal mask - Large	10 pack	
BC803-10	Infant nasal mask - Extra large	10 pack	
BC300-05	Infant bonnet (17-22cm)	5 pack	
BC303-05	Infant bonnet (22-25cm)	5 pack	
BC306-05	Infant bonnet (25-29cm)	5 pack	
BC309-05	Infant bonnet (29-36cm)	5 pack	
BC325-05	Infant head gear (29-36cm)	5 pack	
BC328-05	Infant head gear (35-40cm)	5 pack	
BC331-05	Infant head gear (40-45cm)	5 pack	
BC351-10	Chin strap (20-26cm)	10 pack	
BC353-10	Chin strap (26-32cm)	10 pack	
BC355-10	Chin strap (32-38cm)	10 pack	
BC357-10	Chin strap (38-44cm)	10 pack	










Infant CPAP therapy accessories

Product Code	Description	Quantity	
BC400-10	Gas supply line 2.1m for blender	10 pack	
BC402-05	Gas supply line kit for dual flow meters	5 pack	
BC501-01	Bird® MicroBlender 2-16 L/min DISS	1 each	
BC507-01	Blender oxygen/air Bio-Med 2-16 L/min DISS - White	1 each	
BC510-01	Bird® MicroBlender pole mount bracket	1 each	

Infant nasal high flow breathing circuits

Product Code	Description	Quantity	
RT329	Infant continuous flow circuit (0.3 to 10 L/min) inspiratory heated breathing circuit with autofeed chamber & pressure relief manifold	10 pack	
RT330	Optiflow™ Junior inspiratory heated breathing circuit with autofeed chamber & pressure relief manifold	10 pack	

Infant nasal high flow interfaces

Product Code	Description	Quantity	
OPT312	Optiflow™ Junior nasal cannula - Premature (Wigglepad™ replacement = OPT010)	20 pack	
OPT314	Optiflow™ Junior nasal cannula - Neonatal (Wigglepad™ replacement = OPT012)	20 pack	
OPT316	Optiflow™ Junior nasal cannula - Infant (Wigglepad™ replacement = OPT012)	20 pack	
OPT318	Optiflow™ Junior nasal cannula - Pediatric (Wigglepad™ replacement = OPT012)	20 pack	
BC2425-20	Oxygen therapy nasal cannula - Premature	20 pack	
BC2435-20	Oxygen therapy nasal cannula - Neonatal	20 pack	
BC2745-20	Oxygen therapy nasal cannula - Infant	20 pack	
BC2755-20	Oxygen therapy nasal cannula - Intermediate infant	20 pack	
BC3780-20	Oxygen therapy nasal cannula - Pediatric	20 pack	

Infant nasal high flow accessories


Product Code	Description	Quantity	
OPT010	Optiflow™ Junior Wigglepads™ - Premature (20)	2 pack	
OPT012	Optiflow™ Junior Wigglepads™ - Infant (20)	2 pack	
OPT014	Optiflow™ Junior oxygen tubing	10 pack	
OPT016	Optiflow™ Junior adaptor kit	10 pack	
BC115-05	Pressure relief manifold 40 cm (Blue) suitable for oxygen therapy delivery to infants using systems such as the F&P RT329 or RT330	5 pack	

F&P 810™ System





The F&P MR810 System is designed to deliver airway humidity to match the patient's respiratory and therapy needs. When the MR810 Heated Humidifier is used with the F&P Evatherm™ heated breathing circuit, humidity output performance is enhanced while reducing the unwanted effects of condensate. These benefits are provided with only one probe connection to the circuit, simplifying setup.








Humidifier & accessories

Product Code	Description	Quantity	
MR810XXX	MR810 respiratory heated humidifier (XXX indicates regional plug variation)	1 each	








Evatherm™ breathing circuits

Product Code	Description	Quantity	
900MR810E	Adult Evatherm™ inspiratory heated breathing circuit with expiratory water trap	10 pack	
900MR810	Adult Evatherm™ inspiratory heated breathing circuit	10 pack	



Humidifier brackets & hardware

Product Code	Description	Quantity	
900MR030	Mounting bracket - Pole clamp 15-50mm	1 each	
900MR086	Mounting bracket - Wall	1 each	
900MR087	Mounting bracket - Waterbag pole	1 each	
900MR088	Mounting bracket rail - Dräger/Siemens	1 each	
900MR170	Mounting bracket - C-clamp	1 each	
900MR301	Mounting bracket - Bird 44 mm slot	1 each	
900MR302	Mounting bracket - Concha conversion	1 each	
900MR303	Mounting bracket - Rail	1 each	
900MR304	Mounting bracket - PB840	1 each	
900MR305	Mounting bracket - Dräger Evita 4	1 each	



Humidifier brackets & hardware

Product Code	Description	Quantity	
900MR012	Waterbag pole - 1.0m (40")	1 each	
900MR290	Waterbag pole bent - 0.8 m (28") high	1 each	
900MR291	Waterbag pole bent - 1.08 m (42") high	1 each	
900MR292	Humidifier stand with telescopic waterbag pole - 1.64-2.64m (64-104")	1 each	
900MR293	Humidifier stand with waterbag pole - 1.65m (65")	1 each	
900MR306	Humidifier stand basket	1 each	
900MR599	Oxygen flow adaptor - 22 mm Male to 9/16"	1 each	



Humidifier accessories for reusable systems

Product Code	Description	Quantity	
900MR858	MR850 inspiratory heater wire adaptor	1 each	
900MR859	MR850 inspiratory heater wire adaptor	1 each	










Reusable humidification chambers

Product Code	Description	Quantity	
MR370	Adult humidification chamber	1 each	
MR340X	Infant/Neonatal humidification chamber "X" indicates regional variation	1 each	











Reusable breathing circuit kits

Product Code	Description	Quantity	
900MR784	Adult reusable breathing circuit kit for F&P 850 & 810 Systems	1 each	
900MR782	Neonatal reusable breathing circuit kit for F&P 850 & 810 Systems	1 each	




Adult reusable breathing circuit (900MR784) spare parts

Product Code	Description	Quantity	
900MR126	Adult Y-piece	1 each 50 pack (900MR126L) 100 pack (900MR126C)	
900MR532	Temperature probe housing 22 mm Female - 22 mm Male/15 mm Female	1 each 50 pack (900MR532L) 100 pack (900MR532C)	
900MR534	Tube joint 22 mm Male/15 mm Female	1 each 50 pack (900MR534L) 100 pack (900MR534C)	
900MR715	Reusable adult 1.5 m inspiratory heaterwire with cloverleaf plug	1 each	
900MR139	Autoclavable water trap	1 each 50 pack (900MR139L) 100 pack (900MR139C)	
900MR062	Reusable adult tubing - 0.6m (24")	1 each	
900MR067	Reusable adult tubing - 1.5m (60")	1 each	
900MR074	Reusable adult tubing - 1.1m (44")	1 each	
900MR070	Draw wire - 1.5m (60")	1 each	

Infant reusable breathing circuit (900MR782) spare parts

Product Code	Description	Quantity	
900MR143	Adaptor 10 mm Male to 10 mm Male	10 pack	
900MR145	Neonatal Y-piece pressure port	1 each	
900MR146	Temperature probe housing 10 mm Male-15 mm Male/10 mm Female	1 each 100 pack (900MR146C)	
900MR157	Adaptor 10 mm Male to 15 mm Male/11 mm Female (Drager)	1 each	
900MR178	Adaptor 12 mm Male to 15 mm Male	1 each	
900MR139	Autoclavable water trap	1 each 50 pack (900MR139L) 100 pack (900MR139C)	
900MR406	Adaptor 22 mm Female-10 mm Male	1 each	
900MR711	Reusable infant heaterwire with cloverleaf plug	1 each	
900MR017	Reusable infant tubing - 0.3m (12")	1 each	
900MR025	Reusable neonatal tubing - 0.7 m (28")	1 each	

Infant reusable breathing circuit (900MR782) spare parts

Product Code	Description	Quantity	
900MR026	Reusable neonatal tubing - 1.1m (43")	1 each	
900MR070	Draw wire - 1.5m (60")	1 each	
900MR075	Reusable silicon pressure tubing - 1.65m (66")	1 each	

Index

900IW130X.....	B - 10	900PT421.....	A - 3
900MR012.....	D - 2	900PT422.....	A - 3
900MR017.....	E - 3	900PT426.....	A - 3
900MR025.....	E - 3	900PT427.....	A - 3
900MR026.....	E - 4	900PT501.....	A - 2
900MR030.....	D - 1	900PT531.....	A - 2
900MR062.....	E - 2	900PT600.....	A - 3
900MR067.....	E - 2	900PT601.....	A - 3
900MR070.....	E - 2	900PT602.....	A - 3
900MR070.....	E - 4	900PT603.....	A - 3
900MR074.....	E - 2	900PT913.....	A - 3
900MR075.....	E - 4	900RD008.....	B - 11
900MR086.....	D - 1	900RD009.....	B - 11
900MR087.....	D - 1	900RD010.....	B - 9
900MR088.....	D - 1	900RD014-10.....	B - 9
900MR126.....	E - 2	900RD015-10.....	B - 9
900MR139.....	E - 2	900RD016-10.....	B - 9
900MR139.....	E - 3	900RD101.....	B - 11
900MR143.....	E - 3	900RD110.....	B - 9
900MR145.....	E - 3	900RD301.....	B - 11
900MR146.....	E - 3	900RD302.....	B - 11
900MR148.....	B - 11	BC115-05.....	B - 17
900MR157.....	E - 3	BC151-10.....	B - 12
900MR170.....	D - 1	BC161-10.....	B - 12
900MR178.....	E - 3	BC171-10.....	B - 12
900MR190.....	B - 3	BC190-05.....	B - 13
900MR191X.....	B - 3	BC191-05.....	B - 13
900MR208.....	B - 2	BC192-05.....	B - 13
900MR290.....	D - 2	BC300-05.....	B - 14
900MR291.....	D - 2	BC303-05.....	B - 14
900MR292.....	D - 2	BC306-05.....	B - 14
900MR293.....	D - 2	BC309-05.....	B - 14
900MR301.....	D - 1	BC325-05.....	B - 14
900MR302.....	D - 1	BC328-05.....	B - 14
900MR303.....	D - 1	BC331-05.....	B - 14
900MR304.....	D - 1	BC351-10.....	B - 14
900MR305.....	D - 1	BC353-10.....	B - 14
900MR306.....	D - 2	BC355-10.....	B - 14
900MR406.....	E - 3	BC357-10.....	B - 14
900MR532.....	E - 2	BC400-10.....	B - 15
900MR534.....	E - 2	BC402-05.....	B - 15
900MR599.....	D - 2	BC501-01.....	B - 15
900MR711.....	E - 3	BC507-01.....	B - 15
900MR715.....	E - 2	BC510-01.....	B - 15
900MR782.....	E - 1	BC800-10.....	B - 14
900MR784.....	E - 1	BC801-10.....	B - 14
900MR805.....	B - 2	BC802-10.....	B - 14
900MR806.....	B - 2	BC803-10.....	B - 14
900MR810.....	C - 2	BC2425-20.....	B - 16
900MR810E.....	C - 2	BC2435-20.....	B - 16
900MR817.....	B - 2	BC2745-20.....	B - 16
900MR858.....	E - 1	BC2755-20.....	B - 16
900MR859.....	E - 1	BC3020-10.....	B - 13
900MR860.....	B - 2	BC3520-10.....	B - 13
900MR861.....	B - 2	BC3780-20.....	B - 16
900MR863.....	B - 2	BC4030-10.....	B - 13
900MR868.....	B - 2	BC4540-10.....	B - 13
900MR869.....	B - 2	BC5040-10.....	B - 13
900MR870.....	B - 2	BC5050-10.....	B - 13
900MR961.....	B - 5	BC5550-10.....	B - 13
900PT405.....	A - 3	BC5560-10.....	B - 13
900PT411.....	A - 3	BC6060-10.....	B - 13
		BC6070-10.....	B - 13

Index

BC6570-10.....	B - 13
MR210.....	B - 3
MR210X.....	B - 3
MR225.....	B - 3
MR225X.....	B - 3
MR250.....	B - 3
MR250X.....	B - 3
MR290HFV.....	B - 3
MR290V.....	B - 3
MR290VX.....	B - 3
MR340X.....	E - 1
MR370.....	E - 1
MR810XXX.....	C - 2
MR850XXX.....	B - 2
NC100.....	B - 11
OPT010.....	B - 17
OPT012.....	B - 17
OPT014.....	B - 17
OPT016.....	B - 17
OPT312.....	B - 16
OPT314.....	B - 16
OPT316.....	A - 2
OPT316.....	B - 16
OPT318.....	A - 2
OPT318.....	B - 16
OPT542.....	B - 4
OPT544.....	B - 4
OPT546.....	B - 5
OPT570.....	B - 5
OPT842.....	A - 2
OPT842.....	B - 4
OPT844.....	A - 2
OPT844.....	B - 4
OPT846.....	A - 2
OPT846.....	B - 4
OPT870.....	A - 2
OPT870.....	B - 5
OPT942.....	A - 2
OPT942.....	B - 4
OPT944.....	A - 2
OPT944.....	B - 4
OPT946.....	A - 2
OPT946.....	B - 4
PT101XX.....	A - 2
RD020-01.....	B - 11
RD050-01.....	B - 11
RD803-10.....	B - 10
RD804-10.....	B - 10
RD805-10.....	B - 10
RD806-10.....	B - 10
RD807-10.....	B - 10
RD900XXX.....	B - 10
RD1300-10.....	B - 9
RD1330-10.....	B - 9
RD1340-10.....	B - 9
RD1350-10.....	B - 9
RD1360-10.....	B - 9
RD1370-10.....	B - 9
RT010.....	B - 11
RT013.....	A - 2
RT016.....	B - 7
RT019.....	B - 8
RT020.....	B - 8
RT021.....	B - 8
RT022.....	B - 7
RT023.....	B - 7
RT023.....	B - 8
RT024.....	B - 8
RT028.....	B - 8
RT030.....	B - 11
RT032.....	B - 8
RT040.....	B - 6
RT040L.....	B - 6
RT040M.....	B - 6
RT040S.....	B - 6
RT041L.....	B - 6
RT041M.....	B - 6
RT041S.....	B - 6
RT042.....	B - 6
RT042L.....	B - 6
RT043L.....	B - 6
RT043M.....	B - 7
RT043S.....	B - 7
RT061.....	B - 11
RT062.....	B - 11
RT063.....	B - 11
RT064.....	B - 11
RT065.....	B - 11
RT200.....	B - 7
RT201.....	B - 5
RT202.....	B - 4
RT204.....	B - 7
RT224.....	B - 12
RT225.....	B - 10
RT226.....	B - 10
RT227.....	B - 10
RT228.....	B - 10
RT265.....	B - 10
RT266.....	B - 10
RT267.....	B - 10
RT268.....	B - 10
RT301.....	B - 5
RT319.....	B - 5
RT324.....	B - 12
RT329.....	B - 16
RT330.....	B - 16
RT380.....	B - 7
RT385.....	B - 7
RT913.....	B - 8
RT914.....	B - 8

Disclaimer: Please note that the information in this catalogue (including the product information and images) is summarized and provided for illustrative purposes only. Please refer to the relevant user instructions for more information and confirm details with your local Fisher & Paykel Healthcare representative prior to placing an order. Information is subject to change without notice.

For more information please contact
your local Fisher & Paykel Healthcare representative

Manufacturer 

Fisher & Paykel Healthcare Ltd
15 Maurice Paykel Place
East Tamaki, Auckland 2013

PO Box 14 348 Panmure
Auckland 1741
New Zealand

Tel: +64 9 574 0100
Fax: +64 9 574 0158
Email: info@fphcare.co.nz
Web: www.fphcare.com

Australia (Sponsor)

Fisher & Paykel Healthcare Pty Limited
36-40 New Street, PO Box 167
Ringwood, Melbourne
Victoria 3134, Australia

Tel: +61 3 9879 5022
Fax: +61 3 9879 5232

Austria

Tel: 0800 29 31 23
Fax: 0800 29 31 22

Benelux

Tel: +31 40 216 3555
Fax: +31 40 216 3554

Brazil

Fisher & Paykel do Brasil
Rua Sampaio Viana, 277 cj 21
Paraíso, 04004-000
São Paulo - SP, Brazil

Tel: +55 11 2548 8002

China

代理人/售后服务机构:
费雪派克医疗保健(广州)有限公司,
广州高新技术产业开发区科学城科丰
路31号G12栋301号

电话: +86 20 32053486
传真: +86 20 32052132

Finland

Tel: +358 (0)405 406618
Fax: +46 (0)8 36 6310

France

Tel: +33 1 6446 5201
Fax: +33 1 6446 5221

Germany

Tel: +49 7181 98599 0
Fax: +49 7181 98599 66

India

Tel: +91 80 4284 4000
Fax: +91 80 4123 6044

Irish Republic

Tel: 1800 409 011

Italy

Tel: +39 06 7839 2939
Fax: +39 06 7814 7709

Japan

Tel: +81 3 5117 7110
Fax: +81 3 5117 7115

Korea

Tel: +82 2 6205 6900
Fax: +82 2 6309 6901

Northern Ireland

Tel: 0800 132 189

Russia

Tel and Fax: +7 495 782 21 50

Spain

Tel: +34 902 013 346
Fax: +34 902 013 379

Sweden

Tel: +46 8 564 76 680
Fax: +46 8 36 63 10

Switzerland

Tel: 0800 83 47 63
Fax: 0800 83 47 54

Taiwan

Tel: +886 2 8751 1739
Fax: +886 2 8751 5625

Turkey

Fisher Paykel Sağlık Ürünleri
Ticaret Limited Şirketi,
Alinteri Bulvarı 1161/1 Sokak
No. 12-14, P.O. Box 06371 Ostim,
Ankara, Turkey

Tel: +90 312 354 34 12
Fax: +90 312 354 31 01

UK 

Fisher & Paykel Healthcare Ltd
Unit 16, Cordwallis Park
Clivemont Road, Maidenhead
Berkshire SL6 7BU, UK

Tel: +44 1628 626 136
Fax: +44 1628 626 146

Rest of world

Tel: +64 9 574 0100
Fax: +64 9 574 0158